

DiaTech Oncology

Patient Specific Cancer Testing

Specimen Packing and Shipping Instructions

- Seal the transport tube/container tightly
- Label specimen with institution name, patient name, date and time of collection, and anatomical site of collection
- Place specimen and Ice Packs as shown on the transport box
- Place completed patient information forms and specimen transport box in FedEx plastic “Clinical Pak” envelope (cardboard box for effusion specimen) and seal.
- Complete and affix the FedEx Airbill and 4 copies of the commercial invoice form to the outside of the “Clinical Pak” (cardboard box for effusion specimen). Be sure that the airbill is marked “FedEx International Priority” delivery.
- Call DiaTech Oncology at 1-866-556-5356 or (514) 239-5159 (lab manager) with the FedEx Tracking Number. DiaTech must receive the specimen within 24 hours of collection.
- Specimens sent on Friday must be marked for Saturday delivery.
- Call FedEx to schedule a pickup at 1-800-GOFEDEX (1-800-463-3339).

